



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine storage box for 50 ml tubes, 3 x 3 places, transparent

**€14.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Storage box for 15/50 ml cell culture tubes (Falcon) in bright colors.

Two inserts ensure a firm hold of the tubes in the upper as well as in the lower area.

The boxes are temperature resistant down to  $-80^{\circ}\text{C}$  and have a transparent, milky lid.

5 x 5 places for 15 ml cell culture tubes with max. 17 mm diameter.

3 x 3 places for 50 ml cell culture tubes with max. 30 mm diameter.

## Additional Information

---

No.	2-1923
Manufacturer (Brand)	neoLabLine
EAN	4058072058920
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	stackable/ collapsible with removable louvre insert with closure/lid with louvre insert
Wide	137 mm
Height	128 mm
Length	137 mm
Depth	137 mm
Type closure/ lid	Hinged latch/lid
suitable for	Centrifuge tube
Capacity	9
Squares (WxL)	3 x 3
for vessels with DM up to	30.0 mm
for vessels with volume up to	50 ml
for type (micro) vessel	Centrifuge tube
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

