



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryobox black, 100 places

€16.40
plus VAT &
Shipping

Product Images



Description

Storage box made of polypropylene black, for storage of light-sensitive samples.

100 places for 0.5, 1.5, 2.0 ml reaction tubes and cryotubes.

Material PP black, autoclavable, stackable.

Alphanumeric index on the lid.

Temperature resistant from -90 to +121°C.

For max. 53 mm tube height.

100 openings 12.5 x 12.5 mm

Dimensions: 140 x 140 x 60mm (L x W x H)

Additional Information

No.	7-8040
Manufacturer (Brand)	neoLab
EAN	4058072041953
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with closure/lid with grid insert for different vessel types/volumes for light-sensitive specimens alphanumerically coded
Wide	140 mm
Height	60 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	12.5 x 12.5 mm
Squares (WxL)	10 x 10
for vessels with DM up to	12.5 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

