





neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine cryo storage box for grids, transparent

€5.72
plus VAT &
Shipping



Description

Transparent box made of PP with slip lid for equipping with plastic-coated cardboard grids and rater made of polypropylene for vessels with 12 to 20 mm Ø.

Temperature resistant from -90 to +121°C, autoclavable.

Dimensions: 135 x 135 x 50 mm (L x W x H)

Accessories: PP grid (Order No. 2-2894 ff.)

Additional Information

No.	2-1916
Manufacturer (Brand)	neoLabLine
EAN	4058072058852
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	with closure/lid for different vessel types/volumes
Wide	135 mm
Height	50 mm
Length	135 mm
Depth	135 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
Squares (WxL)	variable
for type (micro) vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

