




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® stainless steel cryobox 136x136x78 mm

**€35.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **Cryomaster® cryobox** made of high-quality **stainless steel** is ideally suited for long-term storage of sensitive samples due to its robustness. Where boxes have to be stored over a longer period of time, frequently frozen or thawed, the stainless steel cryobox is the better alternative to conventional cryoboxes made of cardboard or plastic.

With a basic box size of 136 x 136 mm, the cryoboxes fit into all standard cabinet and chest racks.

- Ideal for cell culture laboratories and culture banks
- Available in three different heights
- For cryotubes with a diameter of 12 mm diameter

## Additional Information

---

No.	4-6021
Manufacturer (Brand)	Cryomaster
EAN	4058072191344
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	silver
Material	Stainless steel
Boxes properties	with special chemical resistance with closure/lid for light-sensitive specimens
Wide	136 mm
Height	78 mm
Length	136 mm
Depth	136 mm
Type closure/ lid	Slip closure
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container

