



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine cryobox coated cardboard, green, 136x136x50mm

€3.90
plus VAT &
Shipping

Product Images



Description

The **green** neoLab cryobox is made of cold-resistant, coated cardboard. Thus, this cardboard box is suitable for sample storage at cryogenic temperatures of up to -86°C . By using different grids, the cryobox can also be used flexibly as a storage box in your laboratory. The slip lid has a practical labeling field on the front side for easy identification of the box or its contents.

Dimensions: 136 x 136 x 50 mm (L x W x H)

Additional Information

No.	2-2697
Manufacturer (Brand)	neoLabLine
EAN	4058072068530
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	green
Material	Cardboard PE coated
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert with labeling field for different vessel types/volumes for light-sensitive specimens
Wide	136 mm
Height	50 mm
Length	136 mm
Depth	136 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
suitable for cooling in	Nitrogen tank
for type (micro) vessel	Reaction vessels Cryogenic tubes

