




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine cryobox coated cardboard, yellow, 136x136x50mm

**€3.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The **yellow** neoLab cryobox is made of cold-resistant, coated cardboard. Thus, this cardboard box is suitable for sample storage at cryogenic temperatures of up to  $-86^{\circ}\text{C}$ . By using different grids, the cryobox can also be used flexibly as a storage box in your laboratory. The slip lid has a practical labeling field on the front for easy identification of the box or the contents.

Dimensions: 136 x 136 x 50 mm (L x W x H)

## Additional Information

---

No.	2-2698
Manufacturer (Brand)	neoLabLine
EAN	4058072068547
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow
Material	Cardboard PE coated
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert with labeling field for different vessel types/volumes for light-sensitive specimens
Wide	136 mm
Height	50 mm
Length	136 mm
Depth	136 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
suitable for cooling in	Nitrogen tank
for type (micro) vessel	Reaction vessels Cryogenic tubes

