



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine cryobox for cell culture tubes coated cardboard, white, 155x155x130 mm

€13.20
plus VAT &
Shipping

Product Images



Description

The neoLab cryobox for cell culture tubes is made of cold-resistant, coated cardboard. Thus, this cardboard box is suitable for sample storage at cryogenic temperatures of up to -86°C . With ration inserts for 15 (6 x 6) and 50 (4 x 4) ml cell culture tubes, this box can be used flexibly as a cryobox or storage box in your laboratory. This cryobox for cell culture tubes is an economical alternative to the conventional, expensive plastic boxes. The slip lid has a practical labeling field on the front side for easy identification of the box or the contents.

Dimensions: 155 x 155 x 130 mm (L x W x H)

Additional Information

No.	2-7095
Manufacturer (Brand)	neoLabLine
EAN	4058072008734
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Cardboard PE coated
Boxes properties	with removable grid insert with closure/lid with labeling field in set for different vessel types/volumes for light-sensitive specimens
Wide	155 mm
Width/inside	148 mm
Height	130 mm
Height inside	128.0 mm
Length	155 mm
Depth	155 mm
Depth inside	148 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	15 ml 50 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
for vessels with DM up to	20.0 mm 30.0 mm
for type (micro) vessel	Centrifuge tube

