



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryobox coated cardboard, white, 136x136x20 mm

€4.17
plus VAT &
Shipping

Product Images



Description

The neoLab cryobox for PCR tubes is made of cold-resistant, coated cardboard. The neoLab cryobox is made of cold-resistant, coated cardboard. Thus, this cardboard box is suitable for sample storage at cryogenic temperatures of up to -86°C. By using different rater inserts, the Kryobox can also be used flexibly as a storage box in your laboratory.

Dimensions: 136 x 136 x 20 mm (L x W x H)

Additional Information

No.	2-2890
Manufacturer (Brand)	neoLabLine
EAN	4058072069322
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Cardboard PE coated
Boxes properties	stackable/pluggable with closure/lid with labeling field for light-sensitive specimens
Wide	136 mm
Height	20 mm
Length	136 mm
Depth	136 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	0.2 ml
suitable for cooling in	Nitrogen tank
Squares (WxL)	16 x 16
for vessels with DM up to	6.0 mm
for type (micro) vessel	PCR reaction vessels

