



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Cryoboxes (PC), 81 places, 53 mm high, Rainbow 4-pack (1 each x blue, green, red, yellow)**

**€58.60  
plus VAT &  
Shipping**

## Product Images



## Description

---

Cryo-box for 9 x 9 cryotubes with 1.2/2.0 ml capacity and 12.5 mm Ø.

The stackable box is in three parts and consists of a white bottom part, a colored compartment insert and a transparent lid with alphanumeric labeling for better identification of the samples.

Material: Polycarbonate

Temperature resistance: -190° to +121°C

Dimensions: 135 x 135 x 53 mm high

**Order No. 7-8049 / Rainbow pack** contains 1 box in each color (blue, red, green, yellow).

## Additional Information

No.	7-8049
Manufacturer (Brand)	neoLab
EAN	4058072042004
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue green red yellow
Material	Polycarbonate (PC)
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	135 mm
Height	53 mm
Length	135 mm
Depth	135 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	2 ml 1.2 ml
suitable for cooling in	Nitrogen tank
Tray size	12.5 x 12.5 mm
Squares (WxL)	9 x 9
for vessels with DM up to	12.5 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

