



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo storage boxes (PC), 81 places, 77 mm high, blue, 6 pcs./pack

**€160.60
plus VAT &
Shipping**

Product Images



Description

Cryobox for 9 x 9 cryotubes with 3 to 4 ml content and 12.5 mm Ø.

The stackable box is in three parts and consists of a white bottom part, a blue compartment insert and a transparent lid with alphanumeric labeling for better identification of the samples.

Material: polycarbonate

Temperature resistance: -190° to +121°C

Dimensions: 135 x 135 x 77 mm high

Additional Information

No.	7-8035
Manufacturer (Brand)	neoLab
EAN	4058072041908
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue transparent (clear) white
Material	Polycarbonate (PC)
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	135 mm
Height	77 mm
Length	135 mm
Depth	135 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	3 ml 4 ml
suitable for cooling in	Nitrogen tank
Tray size	12.5 x 12.5 mm
Squares (WxL)	9 x 9
for vessels with DM up to	12.5 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

