



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo storage box PC, yellow, 9 x 9 places, 96 mm high, 5 pcs./pack

€102.70
plus VAT &
Shipping

Product Images



Description

Cryobox for 9 x 9 cryotubes with 3 to 5 ml content and 12.5 mm Ø.

The stackable box is in three parts and consists of a white bottom part, a yellow compartment insert and a transparent lid with alphanumeric marking for better identification of the samples.

Material: polycarbonate

Temperature resistance: -190° to +121°C

Dimensions: 135 x 135 x 96 mm high

Additional Information

No.	7-8023
Manufacturer (Brand)	neoLab
EAN	4058072041809
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow transparent (clear) white
Material	Polycarbonate (PC)
Boxes properties	stackable/assemblable with removable grid insert with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	135 mm
Height	96 mm
Length	135 mm
Depth	135 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	3 ml 4 ml 5 ml
suitable for cooling in	Nitrogen tank
Tray size	12.5 x 12.5 mm
Squares (WxL)	9 x 9
for vessels with DM up to	12.5 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

