




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene cryobox made of PC, 9 x 9 places for 1.5/2.0 ml

**€73.40** As low  
as  
**€59.52**  
plus VAT &  
Shipping

### Product Images



## Description

---

Cryoboxes made of PC for storage in cryogenic chests and cabinets or in liquid nitrogen.

Temperature resistant from -196°C to +121°C, autoclavable.

With grid division and compartment numbering on the lid.

Compartment size 12 mm, suitable for all common cryotubes with 1.5 or 2.0 ml as well as Eppendorf reaction tubes.

## Additional Information

No.	CD-1525
Manufacturer (Brand)	Nalgene
old neoLab article no.	2-5026
EAN	4058072173500
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	white transparent (clear)
Material	Polycarbonate (PC)
Boxes properties	stackable connectable with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	133 mm
Height	52 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Nitrogen tank
Tray size	12 x 12 mm
Squares (WxL)	9 x 9
for vessels with DM up to	12.0 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes

