



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene cryobox made of PC, 9 x9 places for 5 ml

**€456.96**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Cryoboxes made of PC for storage in cryogenic chests and cabinets or in liquid nitrogen.

Temperature resistant from  $-196^{\circ}\text{C}$  to  $+121^{\circ}\text{C}$ , autoclavable.

With grid division and compartment numbering on the lid.

Compartment size 12 mm, suitable for all common cryotubes with 1.5 or 2.0 ml as well as Eppendorf reaction tubes.

## Additional Information

No.	CD-1527
Manufacturer (Brand)	Nalgene
old neoLab article no.	2-5029
EAN	4058072219390
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	white transparent (clear)
Material	Polycarbonate (PC)
Boxes properties	stackable connectable with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	133 mm
Height	95 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
Tray size	12 x 12 mm
Squares (WxL)	9 x 9
for vessels with DM up to	12.0 mm
for vessels with volume up to	5 ml
for type (micro) vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Nitrogen tank

