



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Cryo storage boxes made of PS, 100 slots, yellow

€58.60
plus VAT &
Shipping

Product Images



Description

Cryobox with 100 places for cryotubes with 1 to 2 ml content and 12.5 mm Ø.

It fits 100 tubes with internal thread or 64 tubes with external thread.

Material: polystyrene, suitable for temperatures from -80° to +80°C, not autoclavable.

Lid transparent with compartment numbering for sample identification.

Dimensions: 133 x 133 x 52 mm high

Additional Information

No.	7-8044
Manufacturer (Brand)	neoLab
EAN	4058072041991
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	yellow transparent (clear)
Material	Polystyrene (PS)
Boxes properties	stackable/connectable with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	133 mm
Height	52 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
Tray size	12.5 x 12.5 mm
Capacity	100 64
Squares (WxL)	10 x 10
for vessels with DM up to	12.5 mm
for vessels with volume up to	1 ml 1.2 ml 1.5 ml 2 ml
for type (micro) vessel	Reaction vessels Cryogenic tubes
suitable for cooling in	Ultra Deep Freezer / Deep Freezer

