



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® cryobox set of 5 with attached lid for
100 cryovials 1.5 ml/ 2.0 ml; green, white, red,
yellow, blue; 133x133x50 mm**

**€32.50
plus VAT &
Shipping**

Product Images



Description

Temperature resistant boxes for use from -196° to $+121^{\circ}\text{C}$. With water-repellent plastic coating. Lid, base and all sides are alphanumerically coded. So you can easily and quickly identify the samples. Hinged lid is attached to the box and prevents accidental misplacement or loss. All boxes are supplied with a removable resistant grid made of PP. After adding ventilation holes and using cardboard grids also suitable for storage in liquid nitrogen.

Additional Information

No.	2-0032
Manufacturer (Brand)	neoLab
EAN	4058072089535
autoclavable at 121°C	No
Color	blue yellow green red white
Material	Cardboard
Boxes properties	with closure/lid with grid insert with labeling field in set for different vessel types/volumes for light-sensitive specimens alphanumerically coded
Wide	133 mm
Height	95 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
Number of parts in set	5
for vessels with volume up to	2 ml
suitable for cooling in	Refrigerator Freezer Ultrat freezer Deep freezer Nitrogen container
Tray size	13 x 13 mm
Squares (WxL)	10 x 10
for vessels with DM up to	12.0 mm

