



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine PCR storage box, pink


Special Price

€6.70 was

€8.45

**plus VAT &
Shipping**

Product Images



Description

Storage rack and cryobox for PCR tubes and strips 0.2 ml.

A special feature of this box are the 14 rows, which make it possible to easily lift in and out even 12-tube strips with lids, as the outer row remains free. The extra flat lid prevents the tubes from falling out.

For strips with 8 or 12 tubes.

14 x 14 rows for max. 196 single vessels or 12 pcs. 8-/12-strips.

Stackable

Alphanumeric labeling on the bottom.

Dimensions: 133 x 133 x 30 mm (W x D x H)

Material: PP, autoclavable up to 121°C

Large quantities on request!

Additional Information

No.	2-1912
Manufacturer (Brand)	neoLabLine
EAN	4058072058814
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky) pink
Material	Polypropylene (PP)
Boxes properties	stackable/assemblable with closure/lid with grid insert alphanumerically coded
Wide	133 mm
Height	30 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	0.2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	7 x 7 mm
Squares (WxL)	14 x 14
for vessels with DM up to	6.0 mm
for type (micro) vessel	PCR reaction vessels

