




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine PCR storage box, green


Special Price

€6.70 was

€8.45

**plus VAT &
Shipping**

Product Images



Description

Storage rack and cryobox for PCR tubes and strips 0.2 ml.

A special feature of this box is the 14 rows, which make it easy to lift even 12-well strips with lids in and out, as the outer row remains free. The extra flat lid prevents the tubes from falling out.

For strips with 8 or 12 tubes.

14 x 14 rows for max. 196 single vessels or 12 pcs. 8-/12-strips.

Stackable

Alphanumeric labeling on the bottom.

Dimensions: 133 x 133 x 30 mm (W x D x H)

Material: PP, autoclavable up to 121°C

Large quantities on request!

Additional Information

| | |
|-------------------------------|--|
| No. | 2-1913 |
| Manufacturer (Brand) | neoLabLine |
| EAN | 4058072058821 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) green |
| Material | Polypropylene (PP) |
| Boxes properties | stackable/assemblable with closure/lid with grid insert alphanumerically coded |
| Wide | 133 mm |
| Height | 30 mm |
| Length | 133 mm |
| Depth | 133 mm |
| Type closure/ lid | Hinged latch/lid |
| suitable for | Micro Vessels / Tubes |
| for vessels with volume up to | 0.2 ml |
| suitable for cooling in | Ultra Deep Freezer/ Deep Freezer |
| Tray size | 7 x 7 mm |
| Squares (WxL) | 14 x 14 |
| for vessels with DM up to | 6.0 mm |
| for type (micro) vessel | PCR reaction vessels |

