



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neoLabLine PCR storage box, blue


Special Price

**€6.70** was

**€8.45**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Storage rack and cryobox for PCR tubes and strips 0.2 ml.

A special feature of this box is the 14 rows, which make it easy to lift even 12-well strips with lids in and out, as the outer row remains free. The extra flat lid prevents the tubes from falling out.

For strips with 8 or 12 tubes.

14 x 14 rows for max. 196 single vessels or 12 pcs. 8-/12-strips.

Stackable

Alphanumeric labeling on the bottom.

Dimensions: 133 x 133 x 30 mm (W x D x H)

Material: PP, autoclavable up to 121°C

Large quantities on request!

## Additional Information

No.	2-1914
Manufacturer (Brand)	neoLabLine
EAN	4058072058838
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky) blue
Material	Polypropylene (PP)
Boxes properties	stackable/assemblable with closure/lid with grid insert alphanumerically coded
Wide	133 mm
Height	30 mm
Length	133 mm
Depth	133 mm
Type closure/ lid	Hinged latch/lid
suitable for	Micro Vessels / Tubes
for vessels with volume up to	0.2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	7 x 7 mm
Squares (WxL)	14 x 14
for vessels with DM up to	6.0 mm
for type (micro) vessel	PCR reaction vessels

