



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Thermo Scientific Table Cooling Container

€434.00
plus VAT &
Shipping

Product Images

Description

Benchtop Cooling Container Protect enzymes, cells, reagents and solutions by maintaining freezing temperatures with the benchtop cooling container, even on the lab bench. Array: 3 x 4 Construction/Design: Non-filled/clear lid Maintained Temp, Hours: Hours,Hours , Hours, <l °C , 5.0 Capacity (metric): 0.2 ml to 0.5 ml Holding Capacity: 12 (16 to 17mm diameter tubes) Temperature Range: 0 °C Type: Benchtop Cooling Container Height (anglo): 7.5 in. Length (anglo): 7.75 in. Width (anglo): 5.5 in. Dimensions (L x W x H): 7.75 x 5.5 x 7.5 in. (197 x 140 x 190 mm) Description: 3x4 compartments for 16- 17mm dia.vials Unit Size: Each

Additional Information

No.	CD-1559
Manufacturer (Brand)	Thermo Scientific
VGKL number	108751164
EAN	4066292478398
Gross weight	1.81
Transport temperature	Room temperature
Country of origin	Mexico
Packaging volume	0.00700000 cc
Customs tariff number	39269097

