



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoBox 50 storage box yellow, PP, 5 x 10 places

**€10.20** As low as  
**€8.92**  
plus VAT &  
Shipping

### Product Images



## Description

---

The neoBox 50 is a practical storage box for the safe storage and transport of 50 reaction tubes or cryovials. Made of durable polypropylene, the neoBox 50 is temperature resistant from  $-90^{\circ}\text{C}$  to  $+120^{\circ}\text{C}$  and autoclavable. The compartments of this storage box are arranged in 5 x 10 rows and alphanumerically labeled.

## Additional Information

No.	2-1941
Manufacturer (Brand)	neoLabLine
EAN	4058072059002
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with closure/lid with louvre insert for different vessel types/volumes
Wide	140 mm
Height	55 mm
Length	86 mm
Depth	86 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
for vessels with volume up to	1.2 ml 2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	12 x 12 mm
Squares (WxL)	5 x 10
for vessels with DM up to	12.0 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes
for type vial/vial	Vial. 12 mm Ø

