



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine neoBox 50 storage box transparent, PP, 5 x 10 places

€10.20 As low as
€8.92
**plus VAT &
Shipping**

Product Images



Description

The neoBox 50 is a practical storage box for the safe storage and transport of 50 reaction tubes or cryovials. Made of durable polypropylene, the neoBox 50 is temperature resistant from -90°C to $+120^{\circ}\text{C}$ and autoclavable. The compartments of this storage box are arranged in 5 x 10 rows and alphanumerically labeled.

Additional Information

No.	2-1966
Manufacturer (Brand)	neoLabLine
EAN	4058072059101
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with closure/lid with louvre insert for different vessel types/volumes
Wide	140 mm
Height	55 mm
Length	86 mm
Depth	86 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
Tray size	12 x 12 mm
Squares (WxL)	5 x 10
for vessels with DM up to	12.0 mm
for vessels with volume up to	1.2 ml 2 ml
for type (micro) vessel	Reaction vessels Cryogenic tubes
for type vial/vial	Vial. 12 mm Ø
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

