



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoBox 81 plus, cryobox, 133x133x53 mm (L x W x H), bottom part red, lid transparent/neutral

**€8.39** As low as  
**€6.04**  
plus VAT &  
Shipping

### Product Images



## Description

---

The **neoBox 81 plus** is the ideal **cryobox** for the safe storage and transport of your valuable samples. With its fixed 9x9 grid, the neoBox 81 plus offers space for 81 cryotubes with a diameter of up to 12.5 mm and a maximum height of 50 mm.

Including perforated plate (2-29239), which is inserted into the lower part. This means that each "sample position" is uniquely numbered and assigned. Thus, all 81 positions in the box are unambiguous and there is no mix-up of specimens.

The lid and the lower part are designed in such a way that the lid is always placed in a clearly oriented position - position lock.

The neoBox 81 plus is made of durable polypropylene and is temperature resistant from -90°C to +121°C and autoclavable.

## Additional Information

---

No.	2-2937
Manufacturer (Brand)	neoLabLine
EAN	4058072168254
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
sterile	No
Boxes properties	with closure/lid with grid insert for different vessel types/volumes alphanumerically coded
Wide	133 mm
Height	53 mm
Length	133 mm
Depth	133 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Squares (WxL)	9 x 9
for vessels with DM up to	13.0 mm

