



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolabLine neoBox 81 plus, cryo boxes made of PP, 133x133x53 mm set of 5 all colors assorted, including perforated plate**

**€36.90** As low  
as  
**€19.20**  
plus VAT &  
Shipping

## Product Images



## Description

---

The **neoBox 81 plus** is the ideal **cryobox** for the safe storage and transport of your valuable samples.

With its fixed 9x9 grid, the neoBox 81 plus offers space for 81 cryotubes with a diameter up to 12.5 mm and a maximum height of 50 mm.

Including perforated plate ( 2-2939), which is inserted into the lower part. Thus each "specimen position" is uniquely numbered and assigned.

Thus all 81 positions in the box are unambiguous and there is no mix-up of specimens.

The lid and the lower part are designed in such a way that the lid is always placed in a clearly oriented position - position lock.

The neoBox 81 plus is made of durable polypropylene and is temperature resistant from -90°C to +121°C and autoclavable.

## Additional Information

No.	2-2938
Manufacturer (Brand)	neoLabLine
EAN	4058072178154
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow green blue red transparent (milky)
Material	Polypropylene (PP)
sterile	No
Boxes properties	with closure/lid with grid insert in set for different vessel types/volumes alphanumerically coded
Wide	133 mm
Height	53 mm
Length	133 mm
Depth	133 mm
Squares (WxL)	9 x 9
for vessels with DM up to	13.0 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer

