



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolabLine neoBox 81 plus, cryo boxes
133x133x53 mm (L x W x H) Set of 5: bottom part
transparent, yellow, green, blue, red; lid
transparent/coded**

€30.70 As low
as
€23.00
plus VAT &
Shipping

Product Images



Description

The **neoBox 81 plus with numerically coded lid** is the ideal **cryobox** for the safe storage and transport of your valuable samples. With its fixed 9x9 grid, the neoBox 81 plus offers space for 81 cryotubes with a diameter of up to 12.5 mm and a maximum height of 50 mm. With the optional perforated plate, which is inserted into the lower part, each "sample position" is clearly numbered and assigned. Thus, all 81 positions in the box are unique and there is no mix-up of specimens.

The lid and the lower part are designed in such a way that the lid is always placed in a clearly oriented position - position lock.

The neoBox 81 plus is made of durable polypropylene and is temperature resistant from -90°C to +121°C and autoclavable.

Additional Information

No.	PK-22948
Manufacturer (Brand)	neoLabLine
EAN	4058072165406
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow green blue red transparent (milky)
Material	Polypropylene (PP)
Boxes properties	with grid insert in set for different vessel types/volumes
Wide	133 mm
Height	53 mm
Length	133 mm
Depth	133 mm
for vessels with volume up to	1.5 ml 2 ml
suitable for cooling in	Refrigerator Freezer Ultrat freezer/ Deep freezer
Squares (WxL)	9 x 9
for vessels with DM up to	13.0 mm

