



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine neoBox 100 storage box yellow, PP, 10 x 10 spaces

€10.40 As low
as
€9.18
plus VAT &
Shipping

Product Images



Description

The neoBox 100 is a practical storage box for the safe storage and transport of 100 reaction tubes or cryovials. The neon-colored base and tight-fitting transparent hinged lid are made of durable polypropylene, making it temperature-resistant from -90°C to $+120^{\circ}\text{C}$ and autoclavable. The compartments of this storage box are arranged in 10 x 10 rows and alphanumerically labeled.

Additional Information

No.	2-1935
Manufacturer (Brand)	neoLabLine
EAN	4058072058951
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with closure/lid with louvre insert for different vessel types/volumes
Wide	140 mm
Height	55 mm
Length	153 mm
Depth	153 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
for vessels with volume up to	1.2 ml 2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Tray size	12 x 12 mm
Squares (WxL)	10 x 10
for vessels with DM up to	12.0 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes
for type vial/vial	Vial. 12 mm Ø

