




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoBox-200 storage box, black

**€35.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Stackable storage box made of PP with 20 x 10 compartments** for tubes up to 12 mm Ø, maximum height 46 mm (HPLC vials, cryotubes, reaction tubes 1.5/2.0 ml).

**Tightly closing hinged lid .**

**Alphanumeric grid on the lid.**

Temperature resistant from -80°C and 121°.

**Material:** Polypropylene, autoclavable up to 121°C

**Dimensions:** 275 x 140 x 60 mm (W x D x H)

Accessories: Insert for reaction vessels 1.5 ml (Order No. 2-1210)

## Additional Information

No.	2-1964
Manufacturer (Brand)	neoLabLine
EAN	4058072059088
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with removable grid insert with closure/lid with grid insert for different vessel types/volumes for light-sensitive specimens alphanumerically coded
Wide	275 mm
Height	60 mm
Length	140 mm
Depth	140 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
for vessels with volume up to	1.2 ml 2 ml
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Squares (WxL)	20 x 10
for vessels with DM up to	12.0 mm
for type (micro) vessel	Reaction vessels Cryogenic tubes
for type vial/vial	Vial. 12 mm Ø

