




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoBox-50, f. Vessels up to 9 mm, 35 mm high, transparent

**€12.10** As low as  
**€10.90**  
plus VAT &  
Shipping

### Product Images



## Description

---

### **Stackable storage box made of PP with 5 x 10 places.**

Temperature resistant from -90°C to 121°C, autoclavable.

For vessels up to 9 mm Ø, max. 35 mm high.

Tightly closing slip lid.

Dimensions: 120 x 64 x 38 mm (W x D x H)

## Additional Information

No.	2-1965
Manufacturer (Brand)	neoLabLine
EAN	4058072059095
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Boxes properties	stackable/connectable with closure/lid with louvre insert for different vessel types/volumes
Wide	120 mm
Height	38 mm
Length	64 mm
Depth	64 mm
Type closure/ lid	Hinged latch/lid
suitable for	Vials Micro Vials/ Tubes
Tray size	9 x 9 mm
Capacity	50
Squares (WxL)	5 x 10
for vessels with DM up to	9.0 mm
for vessels with volume up to	1.2 ml 2 ml
for type (micro) vessel	Reaction vessels Cryogenic tubes
for type vial/vial	Vial. 12 mm Ø
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

