



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® neMultibox red, PP, high form, for all MultiTray tips up to 1,200 µl

€9.00
plus VAT &
Shipping

Product Images



Description

Pipette tip handling made easy!

The practical neoMultibox plus significantly reduces the amount of waste in your laboratory.

To remove the tips, the dust-tight lid of the neoMultibox plus can be opened, pushed open as far as required or removed completely.

The box made of PP can be autoclaved several times at 121°C (20 min.).

Accessories for neoMultiTray pipette tips (order no. 7-3041 to 7-3043)

Additional Information

No.	7-3044
Manufacturer (Brand)	neoLab
EAN	4058072112684
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	red
Material	Polypropylene (PP)
Boxes properties	with closure/lid
Wide	95 mm
Height	80 mm
Length	126 mm
Depth	126 mm
suitable for	Pipette tips (trays and racks)

