



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Beta safety box for inserts, acrylic glass 10 mm thick

**€188.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Beta safety box with hinged lid made of 10 mm thick acrylic glass.

Inner dimensions: 155 x 80 x 40 mm (L x W x H)

## Additional Information

---

No.	4-1011
Manufacturer (Brand)	neoLab
EAN	4058072020507
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polymethylmetacrylate (PMMA) (acrylic glass)
Boxes properties	with removable grid insert with closure/lid with grid insert with beta ray protection for different vessel types/volumes
Wide	100 mm
Height	60 mm
Length	175 mm
Depth	175 mm
Type closure/ lid	Hinged latch/lid
for vessels with volume up to	0.5 ml 1.5 ml 2 ml
suitable for cooling in	Refrigerator
Squares (WxL)	4 x 8
for vessels with DM up to	12.0 mm

