



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Biohazard cooler bag, 32 L

€230.00
plus VAT &
Shipping

Product Images



Description

The **Biohazard cooling bags** are perfectly suitable for the transport of biological material of category B (UN3373, ADR 2013 / P650). The bag wall, reinforced thanks to the polyurethane foam, guarantees an effective isothermal storage of the cooled and heated samples from -20 ° to 37 °C. Likewise, the removable stabilizing wall with a thickness of 3 mm provides optimal cushioning protection against shocks.

- With clear marking - UN3373
- With carrying strap
- Interior with pockets for storing cold packs and documents
- Available in 4 sizes

Additional Information

No.	1-1393
Manufacturer (Brand)	neoLab
EAN	4058072160869
autoclavable at 121°C	No
Color	Blue Grey
Volume	32 l
sterile	No
Boxes properties	insulating
Wide	400 mm
Width/inside	370 mm
Height	360 mm
Height inside	340.0 mm
Depth	400 mm
Depth inside	370 mm
Dimensions	400 x 350 mm

