



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Insulating box, foldable, 23 L

€23.00
plus VAT &
Shipping

Product Images



Description

Practical and lightweight insulated box that can be easily folded after use.

This ensures space savings of up to 60%.

The contents can be stored and transported cold or warm, depending on the application.

Ideally suited for sensitive samples in the laboratory, pharmacy or catering.

Volume: 23 liters

Additional Information

No.	2-1341
Manufacturer (Brand)	neoLab
EAN	4058072087333
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	gray yellow
Material	Expanded polypropylene (EPP)
Volume	23 l
Boxes properties	with closure/lid insulating
Wide	400 mm
Width/inside	360 mm
Form	rectangular
Height	220 mm
Height inside	200.0 mm
Length	400 mm
Depth	400 mm
Depth inside	360 mm
Length inside	360 mm
Type tub/ box	Folding box

