




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Transport container 10 l, 445 x 345 x 90 mm

€24.10
plus VAT &
Shipping

Product Images



Description

Stackable storage and transport container with recessed handles, without side openings.
High chemical and thermal resistance, unbreakable, smooth surface, physiologically harmless, food safe.
Material: Polypropylene (PP)
PP is autoclavable up to 121°C.

Additional Information

No.	2-1301
Manufacturer (Brand)	neoLab
EAN	4058072055950
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Volume	10 l
Boxes properties	stackable/pluggable with special chemical resistance
Wide	345 mm
Width/inside	310 mm
Form	angular
Height	90 mm
Height inside	85.0 mm
Length	445 mm
Depth	445 mm
Depth inside	410 mm
Length inside	410 mm
Type tub/ box	Transport tub/box storage box

