



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Transport container 28 l, 660 x 450 x 130 mm

€46.60
plus VAT &
Shipping

Product Images



Description

Stackable storage and transport container with recessed handles, without side openings.
High chemical and thermal resistance, unbreakable, smooth surface, physiologically harmless, food resistant.
Material: HDPE

Additional Information

No.	2-1303
Manufacturer (Brand)	neoLab
EAN	4058072055974
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	High-density polyethylene (HDPE)
Volume	28 l
Boxes properties	stackable/pluggable with special chemical resistance
Wide	450 mm
Width/inside	400 mm
Form	angular
Height	130 mm
Height inside	120.0 mm
Length	660 mm
Depth	660 mm
Depth inside	600 mm
Length inside	600 mm
Type tub/ box	Transport tub/box storage box

