



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® Safety transport box made of PC, yellow

€172.60
plus VAT &
Shipping

Product Images



Description

Ideal for transporting biological and clinical samples, as well as instruments and products that must be kept clean and dry under test conditions.

The clarity of the neoLab safety transport box allows users to visibly check sample integrity so they can be opened with proper precautions. All boxes are leak-proof, polycarbonate construction. Lid opens to 180 °.

380 x 196 x 160 mm (LxWxH) Outer dimensions

335 x 155 x 135 mm (LxWxH) Inner dimensions

Autoclavable

Additional Information

No.	3-2142
Manufacturer (Brand)	neoLab
EAN	4058072117283
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (yellow)
Material	Polycarbonate (PC)
Wide	196 mm
Height	160 mm
Length	380 mm
Depth	380 mm

