



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® aluminum pipette case 65 x 65 mm, 260/340 mm long

**€41.10** As low as  
**€32.10**  
plus VAT &  
Shipping

### Product Images



## Description

---

Aluminum pipette boxes for storing and sterilizing glass pipettes.  
With slip-over lid.

To protect the pipettes, the base and lid are lined with heat-resistant silicone (up to +160°C).  
Incorporated grooves increase grip.

Material: deep-drawn aluminum  
Basic dimensions: approx. 65 x 65 mm

## Additional Information

---

No.	2-2094
Manufacturer (Brand)	neoLabLine
EAN	4058072059736
Transport temperature	Room temperature
Material	Aluminum
Wide	65 mm
Height	65 mm
Length	260 mm
Depth	260 mm
TBST MAX	160 °C

