



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® aluminum pipette case 65 x 65 mm, 320/470 mm long

€60.30 As low
as
€52.90
plus VAT &
Shipping

Product Images



Description

Aluminum pipette boxes for storing and sterilizing glass pipettes.
With slip-over lid.

The extra-long lid (17 cm) of the pipette box **2-2096** allows an even more variable overall length.
To protect the pipettes, the base and lid are lined with heat-resistant silicone (up to +160°C).
Incorporated grooves increase grip.

Material: deep-drawn aluminum
Basic dimensions: approx. 65 x 65 mm

Additional Information

No.	2-2096
Manufacturer (Brand)	neoLabLine
EAN	4058072059750
Transport temperature	Room temperature
Material	Aluminum
Wide	65 mm
Height	65 mm
Length	320 mm
Depth	320 mm
TBST MAX	160 °C

