



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® stainless steel jar with screw cap, 20 ml

€89.00
plus VAT &
Shipping

Product Images



Description

High-quality small containers made of stainless steel V4A with polished surfaces.

Waterproof even without seal.

Burr-free interior without edges and corners, therefore easy to clean.

Additional Information

| | |
|--------------------------|---|
| No. | 3-2127 |
| Manufacturer (Brand) | neolab |
| EAN | 4058072019136 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| diameter | 38 mm |
| Color | silver |
| Material | Stainless steel |
| Volume | 20 ml |
| sterile | No |
| Form | around |
| Height | 36 mm |
| Type closure/ lid | Screw cap |
| suitable for | variable applications |
| Cans and tins properties | with closure/lid food grade electrically conductive |
| Color closure/ lid | silver |
| Material closure/ lid | Stainless steel |
| Closure / lid properties | closed |
| suitable for cooling in | Refrigerator freezer/ freezer |

