



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoLab® stainless steel jar with screw cap, 100 ml

€189.00
plus VAT &
Shipping

Product Images



Description

High-quality small containers made of stainless steel V4A with polished surfaces.

Waterproof even without seal.

Burr-free interior without edges and corners, therefore easy to clean.

Additional Information

No.	3-2129
Manufacturer (Brand)	neolab
EAN	4058072019150
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	50 mm
Color	silver
Material	Stainless steel
Volume	100 ml
sterile	No
Form	around
Height	73 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with closure/lid food grade electrically conductive
Closure / lid properties	closed
suitable for cooling in	Refrigerator freezer/ freezer
Color closure/ lid	silver
Material closure/ lid	Stainless steel

