



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® screw cap jars 25 ml, 25 pcs./pack

**€17.70**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

White cans made of PP with tightly closing, red screw cap, as sample packaging or for storing powdery or liquid media, small parts, etc.

Autoclavable up to 121°C.

## Additional Information

---

No.	2-2177
Manufacturer (Brand)	neoLab
EAN	4058072060008
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	31 mm
Color	White
Material	Polypropylene (PP)
Volume	25 ml
sterile	No
Form	around
Height	41 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	110 °C
TBST MIN	-10 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	red
Material closure/ lid	Polypropylene (PP)

