



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® screw cap jars 90 ml, 25 pcs./pack**

**€28.90**  
**plus VAT &**  
**Shipping**

## Product Images



## Description

---

White cans made of PP with tightly closing, red screw cap, as sample packaging or for storing powdery or liquid media, small parts, etc.

Autoclavable up to 121°C.

## Additional Information

---

No.	2-2178
Manufacturer (Brand)	neoLab
EAN	4058072060015
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	53 mm
Color	White
Material	Polypropylene (PP)
Volume	90 ml
sterile	No
Form	around
Height	52 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	110 °C
TBST MIN	-10 °C
Color closure/ lid	red
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer

