



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® screw cap jars 1000 ml, 3 pcs./pack

€18.60
plus VAT &
Shipping

Product Images



Description

White cans made of PP with tightly closing, red screw cap, as sample packaging or for storing powdery or liquid media, small parts, etc.

Autoclavable up to 121°C.

Additional Information

No.	2-2182
Manufacturer (Brand)	neoLab
EAN	4058072060053
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	106 mm
Color	White
Material	Polypropylene (PP)
Volume	1000 ml
sterile	No
Form	around
Height	128 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	110 °C
TBST MIN	-10 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	red
Material closure/ lid	Polypropylene (PP)

