



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® screw cap jars made of PS, 47 mm Ø,
70 mm high, 10 pcs./pack**

**€6.11
plus VAT &
Shipping**

Product Images



Description

Polystyrene screw cap jars for clear storage of small parts or powdery media.
Temperature resistance: -40°C to +75°C

Additional Information

No.	2-6983
Manufacturer (Brand)	neoLab
EAN	4058072008376
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	47 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	120 ml
sterile	No
Form	around
Height	70 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	75 °C
TBST MIN	-40 °C
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)

