



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® screw cap jars made of PS, 47 mm  
Ø,110 mm high, 10 pcs./pack**

**€21.80  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Polystyrene screw cap jars for clear storage of small parts or powdery media.

Temperature resistance: -40°C to +75°C

## Additional Information

---

No.	2-6985
Manufacturer (Brand)	neoLab
EAN	4058072008390
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	47 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	190 ml
sterile	No
Form	around
Height	110 mm
Type closure/ lid	Screw cap
suitable for	variable applications
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	75 °C
TBST MIN	-40 °C
Color closure/ lid	transparent (crystal clear)
Material closure/ lid	Polystyrene (PS)
Closure / lid properties	closed
suitable for cooling in	Freezer/ freezer

