



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton glass jar with cap, 500 ml, 12 pcs/pack

€91.80
plus VAT &
Shipping

Product Images



Description

Glass jars made of soda-lime glass type III, particularly resistant.
Ideal for soil samples or retained samples, as the straight wall allows effortless emptying.
The screw cap has a PE insert with Teflon® coating.

Additional Information

No.	WD-1836
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0064
EAN	4058072427948
Transport temperature	Room temperature
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	500 ml
sterile	No
Color bottle	transparent (clear)
Height	95 mm
Type closure/ lid	Screw cap
DM bottle	91 mm
DM Neck	89 mm
Color closure/ lid	white
Thread	89-400 (US)
Neck size	far
Material gasket / septum	Low-density polyethylene (LDPE)
Material bottle	Soda-lime glass (normal glass) 2
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

