



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Nalgene straight wall PPCO wide mouth jars with closure 125 ml Case of 36 70 mm

**€127.08**  
**plus VAT &**  
**Shipping**

### Product Images

---

## Description

Nalgene™ straight wall PPCO wide neck jars with closure Store pastes, solids, samples and laboratory accessories in Thermo Scientific™ Nalgene™ PPCO straight wall wide neck jars. Capacity: 125 ml, 4 oz. Capacity (metric): 125 ml capacity (anglo): 4 oz. Color: Translucent For use with (application): Liquid transfer cap DS2153-0700 Contains: Polypropylene caps Material: PPCO Inside neck diameter (metric): 64 mm Number per box: 36 Number per package: 12 Width (metric): 50 g Outside diameter (Anglo): 2.5 in. Outer diameter (metric): 64 mm Closure size: 70 mm Height (Anglo-American): 2.9 inches Height (metric): 74 mm Unit Size: Case of 36

## Additional Information

No.	CD-0400
Manufacturer (Brand)	Nalgene
VGKL number	228721183
EAN	4066292467927
Gross weight	1.27
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	125 ml
Country of origin	Mexico
Packaging volume	0.02800000 cc
Customs tariff number	39231090

