



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene straight wall PPCO wide mouth jars with closure 250 ml Case of 36 70 mm

€32.16 As low as
€25.98
plus VAT &
Shipping

Product Images

Description

Nalgene™ straight wall PPCO wide neck jars with closure Store pastes, solids, samples, and laboratory supplies in Thermo Scientific™ Nalgene™ PPCO straight wall wide neck jars. Capacity: 250 ml, 8 oz. Capacity (metric): 250 ml Capacity (anglo): 8 oz. Color: Translucent For use with (application): Liquid transfer cap DS2153-0700 Contains: Polypropylene caps Material: PPCO Neck inner diameter (metric): 64 mm Number per carton: 36 Number per package: 6 Type: Straight-walled vessel Outer diameter (Anglo-American): 2.95 in. Outer Diameter (metric): 75 mm Clarity: Translucent Closure Size: 70 mm Closure Type: Cap Height (Anglo): 4.7 in. Height (metric): 119 mm Opening: Wide Shape: Round Unit Size: Case of 36

Additional Information

No.	CD-0401
Manufacturer (Brand)	Nalgene
VGKL number	228721184
EAN	4058072180652
Gross weight	2.498
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	250 ml
Country of origin	Mexico
Packaging volume	0.03400000 ccm
Customs tariff number	39231090

