



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene straight wall PPCO wide mouth jars with closure 30 ml Case of 72 43 mm

€151.92
plus VAT &
Shipping

Product Images

Description

Nalgene™ PPCO Straight Wall Wide Neck Cans with Closure Store pastes, solids, samples and laboratory accessories in Thermo Scientific™ Nalgene™ PPCO Straight Wall Wide Neck Cans. Capacity: 30 ml, 1 oz. Capacity (metric): 30 ml capacity (anglo): 1 oz. Color: Translucent Contains: Polypropylene caps Material: PPCO Neck inner diameter (metric): 33 mm Quantity per case: 72 Quantity per package: 12 Type: Straight-walled vessel Outer diameter (Anglo-American): 1.42 in. Outer Diameter (metric): 36 mm Clarity: translucent Closure Size: 43 mm Closure Type: cap Height (Anglo): 1.9 in. Height (metric): 48 mm Opening: Wide Shape: Round Unit Size: Case of 72

Additional Information

| | |
|-----------------------|--------------------|
| No. | CD-0398 |
| Manufacturer (Brand) | Nalgene |
| VGKL number | 228721181 |
| EAN | 4066292467910 |
| Gross weight | 1.55 |
| Transport temperature | Room temperature |
| Material | Polypropylene (PP) |
| Volume | 30 ml |
| Country of origin | Mexico |
| Packaging volume | 0.01600000 ccm |
| Customs tariff number | 39231090 |

