



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene straight wall PPCO wide mouth jars with closure 500 ml Case of 24 120 mm

€153.12
plus VAT &
Shipping

Product Images

Description

Nalgene™ straight wall PPCO wide neck jars with closure Store pastes, solids, samples, and laboratory supplies in Thermo Scientific™ Nalgene™ PPCO straight wall wide neck jars. Capacity: 500 ml, 16 oz. Capacity (metric): 500 ml capacity (anglo): 16 oz. Color: Translucent Contains: Polypropylene caps Material: PPCO Inside neck diameter (metric): 112 mm Quantity per case: 24 Quantity per package: 6 Type: Straight-walled vessel Outer diameter (Anglo-American): 4.41 in. Outer Diameter (metric): 112 mm Clarity: translucent Closure Size: 120 mm Closure Type: cap Height (Anglo): 3.46 in. Height (metric): 88 mm Opening: Wide Shape: Round Unit Size: Case of 24

Additional Information

No.	CD-0402
Manufacturer (Brand)	Nalgene
VGKL number	228721185
EAN	4066292467934
Gross weight	3.19
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	500 ml
Country of origin	Mexico
Packaging volume	0.03900000 cc
Customs tariff number	39231090

