



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**neolab® shipping container 35 ml, 39 x 55 mm,  
lid enclosed, 500 pcs./pack**

**€122.00  
plus VAT &  
Shipping**



## Description

---

TWIN shipping containers made of selected, particularly clear polypropylene.

With tight screw cap.

With specially rounded bottom for easy removal of even smaller tissue samples.

Stable base.

100 ml can serve as a protective vessel for 50 ml.

## Additional Information

---

No.	2-0130
Manufacturer (Brand)	neoLab
EAN	4058072047511
autoclavable at 121°C	Yes
diameter	39 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	35 ml
sterile	No
Height	55 mm
Type closure/ lid	Screw cap
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	110 °C
TBST MIN	-10 °C
Color closure/ lid	blue
Material closure/ lid	High-density polyethylene (HDPE)

