



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® shipping container 50 ml, 39 x 74 mm,
lid enclosed, 500 pcs./pack**

**€71.90
plus VAT &
Shipping**



Description

TWIN shipping containers made of selected, particularly clear polypropylene.

With tight screw cap.

With specially rounded bottom for easy removal of even smaller tissue samples.

Stable base.

100 ml can serve as a protective vessel for 50 ml.

Additional Information

No.	2-0131
Manufacturer (Brand)	neoLab
EAN	4058072047528
autoclavable at 121°C	Yes
diameter	39 mm
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	50 ml
sterile	No
Height	74 mm
Type closure/ lid	Screw cap
Cans and tins properties	with special chemical resistance with closure/lid
TBST MAX	110 °C
TBST MIN	-10 °C
Color closure/ lid	blue
Material closure/ lid	High-density polyethylene (HDPE)

